BCOMING

D2.3 Antibody detection protocols

Project acronym: BCOMING

Project title: Biodiversity Conservation to Mitigate the risks of emerging infectious diseases

Call: HORIZON-CL6-2021-BIODIV-01

Funded by the European Union

Project no.	101059483
Project acronym:	BCOMING
Project title:	Biodiversity Conservation to Mitigate the risks of emerging infectious diseases
Call:	HORIZON-CL6-2021-BIODIV-01
Start date of project:	01.08.2022
Duration:	48 months
Deliverable title:	D2.3 Antibody detection protocols
Due date of deliverable:	31.01.2023
Actual date of submission:	31.01.2023 then 03.05.2024
Deliverable Lead Partner:	IRD
Dissemination level:	Public

Author list

Name	Organization
Peeters M., Ayouba A, Ramassamy JL.	Institut de Recherche pour le Développement
Cappelle J.	CIRAD
Kaba D., Keita A.	CERFIG
Lagostina L.	HZI
Duong V.	IPC

Document History

Version	Date	Note	Revised by
01	24-07-2023	Draft version	WP leaders
02	03-05-2024	Updated version (clarifications included)	IRD and CIRAD

Disclaimer

The content of the publication herein is the sole responsibility of the publishers and it does not necessarily represent the views expressed by the European Commission or its services.

While BCOMING is funded by the European Union, views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Research Executive Agency (REA). Neither the European Union nor the European Research Executive Agency (REA) can be held responsible for them.

While the information contained in the documents is believed to be accurate, the authors(s) or any other participant in the BCOMING consortium make no warranty of any kind with regard to this material including, but not limited to the implied warranties of merchantability and fitness for a particular purpose.

Neither the BCOMING Consortium nor any of its members, their officers, employees or agents shall be responsible or liable in negligence or otherwise howsoever in respect of any inaccuracy or omission herein.

Without derogating from the generality of the foregoing neither the BCOMING Consortium nor any of its members, their officers, employees or agents shall be liable for any direct or indirect or consequential loss or damage caused by or arising from any information advice or inaccuracy or omission herein.

Executive Summary

Multiplex serological assays

BCOMING is a multidisciplinary project that with partners and teams collecting samples from all terrestrial and aquatic ecosystems, from humans and animals, in Cambodia, Cote d'Ivoire, Guinea and Guadeloupe. It is therefore very important to standardize sampling collection and subsequent laboratory analysis. This deliverable concerns the serological analysis from filoviruses and coronaviruses to test for the presence of antibodies in humans and wildlife.

Multiplex assay to test antibodies against Ebolaviruses

The assay has been adapted and optimized based on the initial version of the assay published by Ayouba and colleagues (1). The current version includes 10 recombinant proteins representing different viral regions (nucleoprotein [NP], 40-kDa viral protein [VP40], and glycoprotein [GP]) from five of the six Ebolavirus species: Zaire Ebolavirus, EBOV; Sudanvirus, SUDV; Bundibugyo, BDBV; Reston, RESTV and Bombali virus, BOMV. The test is standardized to test samples collected as plasma, serum or dried blood spots from humans, bats and rodents. SOPs are available in French and English to be shared with all participating laboratories in BCOMING.

(1). Ayouba A, Touré A, Butel C, Keita AK, Binetruy F, Sow MS, Foulongne V, Delaporte E, Peeters M. Development of a Sensitive and Specific Serological Assay Based on Luminex Technology for Detection of Antibodies to Zaire Ebola Virus. *J Clin Microbiol.* 2016 Dec 28;55(1):165-176. doi: 10.1128/JCM.01979-16. PMID: 27795350; PMCID: PMC5228227.

Multiplex assay to test antibodies against Coronaviruses

The assay has been adapted and optimized based on the initial version of the assay published by Ayouba and colleagues (2). The current version includes 20 recombinant proteins representing spike (SP) and nucleoprotein (NP) from the three beta coronaviruses that are highly pathogenic for humans (SARS-CoV-1, SARS-CoV-2 and MERS-COV), four coronaviruses that cause only mild disease (two beta coronaviruses, OC43, and HKU1; two alpha coronaviruses, 229E and NL63,) and three beta coronaviruses that circulate in bats (HKU9, nobecov subgenus isolated from Rousettus bats in China; HKU3, sarbecov subgenus isolated from Rhinolophus bats in China; 133/2005 merbecov subgenus from Tylonycteris bat in China). The test is standardized to test samples collected as plasma, serum

or dried blood spots from humans, bats and rodents. SOPs are available in French and English to be shared with all participating laboratories in BCOMING.

(2) Ayouba A, Thaurignac G, Morquin D, Tuailon E, Raulino R, Nkuba A, Lacroix A, Vidal N, Foulongne V, Le Moing V, Reynes J, Delaporte E, Peeters M. Multiplex detection and dynamics of IgG antibodies to SARS-CoV2 and the highly pathogenic human coronaviruses SARS-CoV and MERS-CoV. J Clin Virol. 2020 Aug;129:104521. doi: 10.1016/j.jcv.2020.104521.

Transfer of technology and training

The protocols for screening of antibodies to Ebola and Coronaviruses have been transferred to CERFIG and staff from CERFIG has been trained. Additional training is planned in October on coupling of beads.

The protocols for coupling of beads and antibody detection against coronaviruses has been transferred to IPC.

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Annex 4. SOP Luminex_Ebola antibody detection-EN

Annex 5. SOP Luminex_ detection anticorps Coronavirus-FR

Annex 6. SOP Luminex_Coronavirus antibody detection -EN

Explanations on the text highlighted in colors in the Annexes:

The text highlighted in yellow is to indicate that the volumes needed are related to a particular number of plates that can be tested. If the user wants to tests more plates he/she should multiply the volumes highlighted in yellow with the factor needed to prepare the plates needed. The text in blue or red colors are choices of style to draw the attention of the reader to the text.

SOP SC_25e	Luminex® - Couplage des Billes aux antigènes
OBJECTIFS	Instructions pour le couplage de protéines recombinantes aux billes pour la technologie xMAP (Luminex®).
PERSONNEL IMPLIQUÉ	Techniciens Chercheurs
DOCUMENTATION / REFERENCES	<p>Bio-Plex Amine Coupling Kit, Instruction Manual Bio-Rad Laboratories, Inc. 2000 Alfred Nobel Drive, Hercules, CA 94547 USA</p> <p>Development of a Sensitive and Specific Serological Assay Based on Luminex Technology for Detection of Antibodies to Zaire Ebola Virus. Ayouba A et al. Journal of Clinical Microbiology. 2016.</p> <p>Multiplex detection and dynamics of IgG antibodies to SARS-CoV2 and the highly pathogenic human Coronaviruses SARS-CoV and MERS-CoV Ayouba A et al. 2020.</p>
RÉALISÉ PAR	Guillaume Thaurignac Signature : Date :
APPROUVÉ PAR	Martine Peeters et Ahidjo Ayouba Signature : Date :

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SOP SC	Luminex® - bead preparation – Coupling with Antigen
GOALS	Instructions for Coupling Recombinant Proteins to Beads with xMap technology
STAFF INVOLVED	Technicians Researchers
DOCUMENTATION AND REFERENCES	<p>Bio-Plex Amine Coupling Kit, Instruction Manual Bio-Rad Laboratories, Inc. 2000 Alfred Nobel Drive, Hercules, CA 94547 USA</p> <p>Development of a Sensitive and Specific Serological Assay Based on Luminex Technology for Detection of Antibodies to Zaire Ebola Virus. Ayouba A et al. Journal of Clinical Microbiology. 2016.</p> <p>Multiplex detection and dynamics of IgG antibodies to SARS-CoV2 and the highly pathogenic human Coronaviruses SARS-CoV and MERS-CoV Ayouba A et al. 2020.</p>
EDITED BY	Guillaume Thaurignac Signature : Date :
APPROVED BY	Martine Peeters et Ahidjo Ayouba Signature : Date :

I. DEFINITION

Luminex®: *this is a generic term designating a luminometer allowing the demonstration of an interaction between ligands (antigen-antibody, primer-target, hormone-receptor, etc.).*

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2. Materials et equipment

A part of the manipulations is sensitive to light, a dark or light-controlled room is necessary.

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- Recombinant protein (Panel description available in part 4)

For storage and resuspension condition respect supplier's recommendations.

- EDC (réf. Pierce: 22980) (sensitive to light! store at -20°C)
- S-NHS (réf. Pierce : 24510) (sensitive to light ! store at -20°C)
- Vortex
- Ultrasonic bath
- Pipets (P10, P100-P200, P1000)
- Magnetic rack (surebeads magnetic rack 16 tube holder)
- rotor
- Eppendorf 1,5mL « Protein LoBind »
- Amine Coupling Kit, Bio Rad
- MagPlex C microspheres (different region), Bio-Rad/Luminex (sensitive to light! store at +4°C) (réf. MagPlex™ MC100XX-01)

3. Manipulation

3. 1. Coupling of recombinant proteins to beads **(To do 6 plates)**

3. 1. 1 Beads *Activation*

- a. Select the colors to activate; identify the tubes by the name of the recombinant protein to be coupled and the bead number
- b. Vortexer beads during 30 sec, sonicate during 30 sec.
- c. Pipet **100 μ l*** of beads in a 1.5mL Eppendorf tube of "Protein LoBind" quality which prevents the beads from sticking to the walls and thus promotes the formation of a clear pellet. By default, use a tube provided in the coupling kit but be very careful about the dispersion of the beads and their not very hermetic closure system (loss during vortexing, evaporation during storage).
- d. Put on the magnetic rack during 1 min.
- e. Gently remove the supernatant (with a p100)
- f. Add 200 μ l of « **bead wash buffer** » ; Vortex 30 sec et sonicate 10 sec ;
- g. Put on the magnetic rack during 1 min.
- h. remove the supernatant (with a p100)
- i. Suspend beads with **80 μ l*** of « **bead activation buffer** » ; Vortex 30 sec and sonicate 30 sec.
- j. Keeping the vortex in rotation, add **10 μ l*** of EDC at 50mg/ml, Vortex immediately afterwards, then very quickly add **10 μ l*** of S-NHS at 50mg/ml. Next Vortex 10 sec.
- k. Vortex 30 sec and incubate (rotor), at 25 rotate by minute during 20min at room temperature. Protect from light during this step with aluminum
- l. After incubation add 150 μ l of PBS (provide with coupling kit)
- m. Vortex 10 sec.
- n. Put on the magnetic rack during 1 min.
- o. remove the supernatant then add 500 μ l of PBS. Vortex 10 sec and put on the magnetic rack during 1 min.
- p. Gently remove the supernatant and et suspend in 100 μ l of PBS Vortex gently during 30 sec and sonicate 15 sec.

3. 1. 2 *Antigen coupling*

- a) Add optimal quantity of antigen (cf panel tab in part 4) to the beads
Adjust volume to 400 µl qsp with PBS 1X.
- b) Incubate during 2h at room temperature at 25 rotate by minute during 2h at room temperature. Protect from light during this step with aluminum (rotor)
- c) Put on the magnetic rack during 1 min ; Gently remove the supernatant (with a p1000)
- d) Wash with 500 µl of PBS ; Vortex 10 sec Put on the magnetic rack during 1 min ; Gently remove the supernatant (with a p1000).
No sonicate.

3. 1. 3 *Saturation*

- a) Suspend in 250µl of « **blocking buffer** » Vortex carefully during 15 sec.
- b) Incubate (rotor), at 25 rotate by minute during 30min at room temperature. Protect from light durin this step with aluminium.
- c) Put on the magnetic rack during 1 min ; Gently remove the supernatant (with a p1000).
- d) Add 500 µl of « **Storage buffer** », Vortex 10 sec then put on the magnetic rack during 2 min ; Gently remove the supernatant (with a p1000).
- e) Suspend beads in **150 µl*** of « **Storage buffer** ».

NB: it is possible to provide an additional "dead" volume (+15µL maximum) which will not alter the reactions and which makes it possible to anticipate volume losses caused by repeated pipetting. Date each tube of paired beads. According to Biorad, they can be stored at 4°C away from light for 1 year. However, it is advisable to validate their performance 3 months after coupling.

The quantities of reagents given for the coupling correspond to a coupling to perform 6 plates. You can adapt those quantities to the real numbers of plate you need. In this case, increase the volumes highlighted in yellow proportionally. The maximum volume to start (volume of native Beads) is 500 (coupling to perform 30 plates)

3. 2. *Coupling Validation*

Before any use of the coupling, it must be validated on a panel of negative and positive controls, respecting the conditions of use described in the SOPs describing the various Screenings.

This validation is essential because it alone guarantees the effectiveness of the manipulations that will be carried out afterwards.

4. Antigen reference and quantity needed

For storage and resuspension condition respect supplier's recommendation.

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Multiplex Immunoassay Luminex® Ebolavirus	
GOALS	Instructions for performing the multiplex serological test using xMAP MagPix® (Luminex®) technology
STAFF INVOLVED	Technicians Researchers
DOCUMENTATION AND REFERENCES	Ayouba A, Touré A, Butel C, Keita AK, Binetruy F, Sow MS, Foulongne V, Delaporte E, Peeters M. Development of a Sensitive and Specific Serological Assay Based on Luminex Technology for Detection of Antibodies to Zaire Ebola Virus. J Clin Microbiol. 2016 Dec 28;55(1):165-176. doi: 10.1128/JCM.01979-16. PMID: 27795350; PMCID: PMC5228227.
EDITED BY	Guillaume Thaurignac and Antoine Nkuba Signature : Date :
APPROVED BY	Martine Peeters et Ahidjo Ayouba Signature : Date :

I. DEFINITION

Luminex®: this is a generic term designating a luminometer allowing the demonstration of an interaction between ligands (antigen-antibody, primer-target, hormone-receptor, etc.).

BSA : Bovine Serum Albumin

PBS : Phosphate Buffered Saline

FBS : foetal bovine sera

II. PRINCIPLE

Luminex® technology works according to the combined principles of ELISA and flow cytometry through the use of fluorescent beads. With the MagPix up to 50 interactions can be studied simultaneously. The beads each have a different color to identify the antigen to which they have been coupled. It is a biotinylated secondary antibody (directed against human IgG for example) which, combined with streptavidin coupled to Phycoerythrin (SAPE), makes it possible to quantify the presence or absence of interaction by an MFI value (Median Fluorescence Intensity).

III. INSTRUCTIONS AND PROCEDURES

1. Security Procedures

All biological material should be considered potentially infectious. Therefore, wearing gloves and a lab coat is mandatory during handling.

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A part of the manipulations are sensitive to light, a dark or light-controlled room is necessary.

Materials et reactive

- A complete **xMap** system : The **MagPix** and it computer with Xponent instaled
- Manual magnetic washer for 96 well plate
- Tablet of 1X PBS
- Tween-20 (polyoxysorbitane monolaurate)
- BSA (SIGMA : A7906)
- Foetal bovine sera (inactivated 30 min 56°C), Gibco
- Biotin Mouse anti-human IgG (BD pharmigen réf. 555785)
- *Biotin goat anti-bat IgG H/L (Bethyl A140-118B)*
- *Biotin goat anti-mouse IgG Y-chain specific (Sigma B7022)*
- *Biotin goat anti-rat IgG (whole molecule) (Sigma B7139)*
- Streptavidin-R-phycoerythrin (SAPE) (Fisher Scientific/Life technologies)
- NaH₂PO₄ (powder)
- Na₂HPO₄ (powder)
- Na Cl (powder)
- 96 well plate with flat bottom, Greiner (réf. 655906) with lid
- Reagent reservoirs
- Vortex

- Plate Agitator
- Pipets (P10, P100-P200, P1000) and without filter compatible tips
- Multichannel Pipets 30 – 300 µL (8 and/or 12 Channel)
- Tube 5mL for mixes preparation
- Tube 1,4mL (micronic) for dilutions with compatible racks
- MagPix calibration kit
- MagPix performance verification kit
- MagPix Drive fluid
- Magnetic bead coupled with antigen
(Avoid expose them to light)
- Isopropanol 70%
- Soude 0.1N
- Bleash
- Refrigerator and freezer
- Magnetic stirrer and bar magnet
- Glass bottle 0.5L, 1L et 2L
- A source of clean water is necessary (distillator, osmosis system)
- Personal protective equipment (gloves, lab coat)
- Ethanol 70% (for cleaning surface)
- heating block or / and laboratory oven
- Microtube 1.5mL and 2mL (and compatible cryoboxe)

List of Antigen

Ag name	suplier	reference	Quantity for size 1X
NP EBOV	Sino Biological	40443-V07E1	2µg
GP EBOV Kiss	Sino Biological	40442-V08B1	2µg
GP EBOV May	Sino Biological	40304-V08B1	2µg
VP40 EBOV	Sino Biological	40446-V07E	2µg
NP SUDV	Sino Biological	40444-V07E	2µg
GP1 SUDV	Sino Biological	40094-V08B	1µg
VP40 SUDV	Sino Biological	40447-V07E	2µg
GP BDBV	Sino Biological	40368-V08B	2µg
VP40 BDBV	Sino Biological	40448-V07E	2µg
GP RESTV	IBT Bioservice	0504-015	2µg
GP1 BOMV	Sino Biological	IT1-2	2µg

3. Procedures

3.1 Buffer preparation

Record the product name and preparation date of each buffer.

- Hypertonic PBS: example for 1 L
NaH₂PO₄= 0,8g
Na₂HPO₄= 2,5g
Na Cl= 88g
Qsp 1L DIH₂O
Store at +4°C.
- 1X PBS: Put 1 tablet in DIH₂O. (respect supplier recommendation)
Store at +4°C.
- Reading buffer: PBS 1X + 1% BSA
Store at +4°C.
- Dilution buffer: (Store at +4°C):
50% hypertonic PBS
1% BSA
5% inactivated FBS (incubate 56°C for 30min)
0.2% Tween-20
Qsp DIH₂O

3.2 Process

3.2.1 Sample preparation

The samples to be tested with Luminex must not present any infectious risk. Before any reading, they must be inactivated in the BSL3 laboratory according to the required health and safety rules.

- if necessary, aliquot in a biosafety atmosphere.

For plasma, 50µL is a good volume for performing several multiplex immunoassays.

- Inactivate **30min to 56°C in heating block** or 1H in a 56°C laboratory oven

After inactivation store aliquot to -20°C

3.2.1 Multiplex Immunoassay with Luminex MagPix

This procedure includes an overnight incubation and therefore takes place over two days.

3.2.1.1 Day 1

1 Sample dilution

Dilutions are carried out in 1.4mL tubes. First of all, it is therefore necessary to prepare a rack with as many 1.4mL tubes as samples to be tested. During this step, the plan of the plate should be produced, which remains identical until the end of the manipulation.

For any reading, remember to include positive controls and controls that will cover the recombinant protein panels used as much as possible and to provide a blank well (buffers alone).

a) Dilution of human plasma and sera

Human plasma/sera samples are tested diluted **1:1000** in sample dilution buffer.

b) Dilution of the eluates (or recovery) of human DBS

On the same principle as above, dilute the DBS eluates to 1/1000 in dilution buffer. The DBS eluates are already diluted to 1/40 after elution so only a 1/25 dilution is necessary

NB - We advise depositing 5µL of eluate in 120µL of dilution buffer.

c) Dilution of bat and rodent plasma and sera

*bat plasma/sera samples are tested diluted 1:2000 in sample dilution buffer.
rodent plasma/sera samples are tested diluted 1:200 in sample dilution buffer.*

d) Dilution of the eluates (or recovery) of bat and rodent DBS

On the same principle as above, dilute the bat's DBS eluates to 1/2000 in dilution buffer. The bat DBS eluates are already diluted to 1/50 after elution so only a 1/40 dilution is necessary

NB - We advise depositing 5µL of eluate in 195µL of dilution buffer.

For rodent dilute the DBS eluates to 1/200 in dilution buffer. The rodent DBS eluates are already diluted to 1/50 after elution so only a 1/4 dilution is necessary

NB - We advise depositing 30µL of eluate in 90µL of dilution buffer.

2 Beads mix preparation

For each of the antigens (and therefore each color of beads (# Number), between 5000 and 8000 beads are required per well. This corresponds to 0.25 µL of each paired bead per well, all in a final volume of 50 µL per well.

For example, if 11 different recombinant proteins have been successfully coupled, the volume of the bead mix for one well will be $11 \times 0.25 = 2.75 \mu\text{L}$ for one well, which will be completed with $47.25 \mu\text{L}$ of dilution buffer.

Then calculate the volume to be pipetted for the number of wells to be made, including a margin.

For example, for a plate (96 wells) it is necessary to prepare a mix for 100 wells.

For 100 wells, $25 \mu\text{L}$ of each paired bead must be deposited per well, all in a final volume of $5000 \mu\text{L}$ qsp of dilution buffer.

Thus, if 11 different recombinant proteins have been successfully coupled, the volume of the mix of beads for a 96-well plate will be $11 \times 25 \mu\text{L} = 275 \mu\text{L}$ which will be completed with $4725 \mu\text{L}$ of dilution buffer.

To prepare the mix of beads, Vortex the vials of coupled beads for 30 sec. Take the required volume of each bead and place them in the sample dilution buffer then Vortex for 30 sec.

Wrap the mix tube in aluminum to protect it from light. The mix, if prepared in advance, can be kept for a maximum of 2 hours at 4°C away from light.

3 Deposit of beads mix and sample dilution

Beforehand, place the plate shaker at 4°C (cold rooms / refrigerator).

Once the sample dilutions and bead mix are ready, take a 96 plate for Luminex® with flat bottom wells.

- Identify and mark the plate.
- Vortex the Mix of Marbles for 30 sec and pour it into a tank (provided for this purpose only)
- Add $50\mu\text{L}$ of bead mix/well to all wells with a multichannel pipette

- Place the plate on the magnetic support for 2 min and wash by inversion

- Distribute $100 \mu\text{L}$ /well of diluted samples by mixing by aspiration/discharge with the multichannel pipette or $100 \mu\text{L}$ of buffer for the blank well.

- Replace the plastic lid on the plate and cover it with aluminum foil. Place it on the stirrer which has been placed at 4°C .
- Set the agitator to 400rpm.

- Leave to incubate for at least 16 hours (overnight) at 4°C (cold room).

3.2.1.2 Day 2

The next morning take out the buffer at room temperature.

1 Getting Started the MagPix

- Put the switch on the back of the device to ON
- Switch on the computer and the MagPix (button on the front)
- Launch the Xponent application (shortcut on the desktop) and identify yourself.
- In the AutoMaint tab select Fluidics Prep
- Fill the tank with 70% Isopropanol in the locations indicated on the screen by the AF zone.
- Click on "Start" to launch the "Fluidics Prep" which lasts approximately 2 minutes.

2 Washing and depositing the anti IgG conjugate

Biotin anti-human IgG-gamma chain specific (500µg/mL) is used at **4µg/mL** (for a final volume of 50µL of conjugate solution/well).

For bat, the secondary antibody is used at 0.1µg/mL (for a final volume of 50µL of conjugate solution/well).

For rodent use a mix of two secondary antibody (anti Rat and anti-Mouse). Mix rat and mouse secondary antibody as a mix. Twice have a final dilution fixed at 1/4000(for a final volume of 50µL of conjugate solution/well).

- Prepare the final conjugate solution in dilution buffer.
- Remove the stirrer plate and remove the aluminum (place the stirrer at room temperature)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 150µL/well of dilution buffer (Washing 1/5)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (2/5 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (3/5 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (4/5 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (5/5 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Off the support, deposit 50µL/well of the final conjugate solution
- Replace the plastic lid on the plate and place it on the stirrer and cover it with aluminum foil.
- Holding the plate firmly, set the stirrer at 1200 RPM for 1 minute to resuspend the beads.
- Reduce the stirrer speed to 400 RPM and incubate for **30 minutes at room temperature.**

3 Performance verification, Calibration and Batch creation

These steps are to be carried out during the 30-minute incubation.

- **MagPix Calibration should be performed weekly when the MagPix is in use.**
 - In the AutoMaint tab select Calibration Verification
 - Fill the tank with 70% Isopropanol in the locations indicated on the screen by the AF zone.
 - Vortex for 30s the 3 vials of beads from the verification kit and the vial from the calibration kit
 - Add 5 to 6 drops of beads of each color in the wells of the bar indicated on the screen
 - Open the hatch from the software using the "eject/retract" button
 - Place the tank and the bar and click again on "eject/retract" to close the hatch.
- Click on "Start" to launch the "Calibration Verification" which lasts about 12 minutes.

NB - This procedure includes checking the performance of the MagPix.

- **The MagPix Performance verification should be performed daily when the Magpix is in use**
 - In the AutoMaint tab select the Performance Verification
 - Fill the tank with 70% Isopropanol in the location indicated on the screen by the AF zone.
 - Vortex for 30s the 3 vials of beads from the verification kit
 - Add 5 to 6 drops of beads of each color in the wells of the bar indicated on the screen
 - Open the hatch from the software using the "eject/retract" button
 - Place the tank and the bar and click again on "eject/retract" to close the hatch.
 - Click on "Start" to launch the "Performance Verification" which lasts about 5 minutes.
- **batch preparation (with an existing protocol)**
 - in the batches tab click on "new batch from an existing protocol" and select "PanCoV protocol".
 - Set the plate format (position of blank and samples)
 - Import the plate plan (linear .txt format) and complete the dilution column
 - In the "post batch routine" pane, select the "post batch routine" command.
 - Fill the tank with distilled water, 70% Isopropanol and 0.1N soda (NaOH) in the locations indicated on the screen.
 - Open the hatch from the software using the "eject/retract" button
 - Place the tank and the bar and click again on "eject/retract" to close the hatch.
 - Save the batch which then appears as pending (Pending)

2 Washing and depositing the Streptavidin-PE (SAPE)

Streptavidin-PE (1mg/mL) is used at **4µg/mL final** (for a final volume of 50µL of solution/well).

Example for a plate: you need 100 wells x 50µL of solution at 4µg/ml = 5 ml of conjugate

solution at 4µg/ml final (i.e. 20 µL of Streptavidin-PE at 1mg/ml) in 4980µL of dilution buffer)

After the incubation with the conjugate:

- Prepare the Streptavidin-PE solution at 4 µg/mL final in dilution buffer
- Cover the tube of diluted Streptavidin-PE with aluminum to protect it from light.
- Remove the agitator plate and remove the aluminum
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 150µL/well of dilution buffer (Washing 1/5)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (2/5 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (3/5 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (4/5 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (5/5 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Off the support, deposit 50µL/well of the Streptavidin-PE solution at 4 µg/mL final
- Replace the plastic lid on the plate and place it on the stirrer and cover it with aluminum foil.
- Holding the plate firmly, set the stirrer at 1200 RPM for 1 minute to resuspend the beads.
- Reduce the stirrer speed to 400 RPM and incubate for 10 minutes at room temperature.

3 Finals wash and plate preparation for reading

After the incubation with the Streptavidin-PE:

- Remove the agitator plate and remove the aluminum
- Place the plate on the magnetic support for 2 min and wash by inversion
- Out of support, deposit 150µL/well of **reading buffer (Washing 1/5)**
- Place the plate on the magnetic support for 2 min and wash by inversion
- Out of support, deposit 100µL/well of **reading buffer (Washing 2/5)**
- Place the plate on the magnetic support for 2 min and wash by inversion
- Out of support, deposit 100µL/well of **reading buffer (3/5 wash)**
- Place the plate on the magnetic support for 2 min and wash by inversion
- Place the plate on the magnetic support for 2 min and wash by inversion
- Out of support, deposit 100µL/well of **reading buffer (4/5 wash)**
- Place the plate on the magnetic support for 2 min and wash by inversion
- Place the plate on the magnetic support for 2 min and wash by inversion
- Out of support, deposit 100µL/well of **reading buffer (5/5 wash)**
- Place the plate on the magnetic support for 2 min and wash by inversion
- Excluding support, deposit 50µL/well of reading buffer
- Replace the plastic lid on the plate and place it on the stirrer and cover it with aluminum foil.
- Holding the plate firmly, set the stirrer at 1200 RPM for 1 minute to resuspend the beads.
- Add 50µl/well of reading buffer to obtain a final volume of 100µl/well.
- Open the hatch from the software using the "eject/retract" button
- Remove the cover and place the plate on the MagPix support and click again on

"eject/retract" to close the hatch

- Open the batch selecting it and clicking on Edit
- Click on the "START" button, and validate so that the reading of the plate starts

Reading a 96-well plate can take approximately 60 minutes. The Post batch routine is triggered automatically at the end of the plate.

4 Result exportation and validation

At the end of the reading, once the post batch routine procedure is finished in the "result" tab, export the results in CSV format (.csv).

This file should be kept as is and should be saved in a results folder on the MagPix computer.

Use a USB flash drive to transfer a copy of this file to another computer for analysis.

For the validation of the plate, several indicators must be checked:

- The consistency of the values obtained with the positive and negative controls with their reference values
- The number of unique beads analyzed for each color of beads (each antigen) in each well (Count). An error is raised if a count is less than 50. The results of a sample that has at least a count less than 25 will not be usable and the sample will have to be ironed.

After validation of the plate, it is possible to integrate the results obtained in Net MFI into a global results file (database)

5 MagPix shut down

In the AutoMaint tab select "shut down"

- Fill the tank with mQ water and bleach diluted at 10% in the places indicated on the screen.
 - Open the hatch from the software using the "eject/retract" button
 - Place the tank and the bar and click again on "eject/retract" to close the hatch.
 - Click on "Start" to launch the shutdown which lasts about 5 minutes.
-
- Turn off the device via the software then de-identify
 - Close the Xponent software by clicking on Exit
 - Set the switch on the back of the device to off
 - Turn off the computer

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Protéines recombinantes	Fournisseurs	références	Quantités pour un couplage pour 6 plaques
NCP-CoV(2019-nCoV) nucleocapsid prot (His tag)	sino Biological	40588-V08B	2µg
SARS-CoV-2(2019-nCoV) Spike protein (RDFB, His tag)	sino Biological	40592-V08B	1µg
SARS-CoV Nucleoprotéin /NP protein (His Tag)	sino Biological	40143-V08B	2µg
SARS-CoV Spike protein (RDFB, His tag)	sino Biological	40150-V08B2	1µg
MERS-CoV Nucleoprotéin /NP protein (His Tag)	sino Biological	40068-V08B	2µg
MERS-CoV Spike/S1 protein (His Tag)	sino Biological	40069-V08B1	1µg
Human coronavirus (HCoV-229E) Spike/S1 protein (S1 Subunit His Tag)	sino Biological	40601-V08H	1µg
Human coronavirus HKU1 (isolate N5) (HCoV-HKU1) Spike/S1 protein (S1 Subunit His Tag)	sino Biological	40602-V08H	1µg
Human coronavirus (HCoV-OC43) Spike protein (S1ECD,His Tag)	sino Biological	4067-V08H1	1µg
Human coronavirus (HCoV-NL63) Spike protein (S1 Subunit His Tag)	sino Biological	40600-V08H	1µg
Recombinant human coronavirus 229E Nucleoprotein(N)	Cusabio	CSB-EP322753HIT	2µg
Recombinant human coronavirus HKU1 Nucleoprotein(N)	Cusabio	CSB-EP607323HIW	2µg
Recombinant human coronavirus OC43 Nucleoprotein(N)	Cusabio	CSB-EP333798HIW	2µg
Recombinant human coronavirus NL63 Nucleoprotein(N)	Cusabio	CSB-EP744150HIX	2µg
Recombinant Bat coronavirus 133/2005 Nucleoprotein(N)	Cusabio	CSB-BP605318BFA	2µg
Recombinant Bat coronavirus 133/2005 Spikevglycoprotein (S) ,partial	Cusabio	CSB-BP610383BFA	1µg
Recombinant Bat coronavirus HKU3 Nucleoprotein(N)	Cusabio	CSB-BP664686BFD	2µg
Recombinant Bat coronavirus HKU3 Spikevglycoprotein (S) ,partial	Cusabio	CSB-BP663395BFD	1µg
Recombinant Bat coronavirus HKU9 Nucleoprotein(N)	Cusabio	CSB-BP385259BOV	2µg
Recombinant Bat coronavirus HKU9 Spikevglycoprotein (S) ,partial	Cusabio	CSB-BP385258BOV	1µg

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Multiplex Immunoassay Luminex® Coronavirus	
GOALS	Instructions for performing the multiplex serological test using xMAP MagPix® (Luminex®) technology
STAFF INVOLVED	Technicians Researchers
DOCUMENTATION AND REFERENCES	Multiplex detection and dynamics of IgG antibodies to SARS-CoV2 and the highly pathogenic human Coronaviruses SARS-CoV and MERS-CoV Ayoub A et al. 2020.
EDITED BY	Guillaume Thaurignac and Antoine Nkuba Signature : Date :
APPROVED BY	Martine Peeters et Ahidjo Ayoub Signature : Date :

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Materials et reactive

- A complete **xMap** system : The **MagPix** and it computer with Xponent installed
- Manual magnetic washer for 96 well plate
- Tablet of 1X PBS
- Tween-20 (polyoxysorbitane monolaurate)
- BSA (SIGMA : A7906)
- Foetal bovine sera (inactivated 30 min 56°C), Gibco
- Biotin Mouse anti-human IgG (BD pharmigen réf. 555785)
- [Biotin goat anti-bat IgG H/L \(Bethyl A140-118B\)](#)
- [Biotin goat anti-mouse IgG Y-chain specific \(Sigma B7022\)](#)
- [Biotin goat anti-rat IgG \(whole molecule\) \(Sigma B7139\)](#)
- Streptavidin-R-phycoerythrin (SAPE) (Fisher Scientific/Life technologies)
- NaH₂PO₄ (powder)
- Na₂HPO₄ (powder)
- Na Cl (powder)
- 96 well plate with flat bottom, Greiner (réf. 655906) with lid
- Reagent reservoirs
- Vortex

- Plate Agitator
- Pipets (P10, P100-P200, P1000) and without filter compatible tips
- Multichannel Pipets 30 – 300 μ L (8 and/or 12 Channel)
- Tube 5mL for mixes preparation
- Tube 1,4mL (micronic) for dilutions with compatible racks
- MagPix calibration kit
- MagPix performance verification kit
- MagPix Drive fluid
- Magnetic bead coupled with antigen
(Avoid expose them to light)
- Isopropanol 70%
- Soude 0.1N
- Bleash
- Refrigerator and freezer
- Magnetic stirrer and bar magnet
- Glass bottle 0.5L, 1L et 2L
- A source of clean water is necessary (distillator, osmosis system)
- Personal protective equipment (gloves, lab coat)
- Ethanol 70% (for cleaning surface)
- heating block or / and laboratory oven
- Microtube 1.5mL and 2mL (and compatible cryoboxe)

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a) Dilution of human plasma and sera

Human plasma/sera samples are tested diluted **1:200** in sample dilution buffer.

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On the same principle as above, dilute the DBS eluates to **1/200** in dilution buffer. The DBS eluates are already diluted to 1/40 after elution so only a 1/5 dilution is necessary

NB - We advise depositing 25µL of eluate in 100µL of dilution buffer.

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*bat plasma/sera samples are tested diluted **1:400** in sample dilution buffer.
rodent plasma/sera samples are tested diluted **1:400** in sample dilution buffer.*

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On the same principle as above, dilute the bat's DBS eluates to 1/400 in dilution buffer. The bat DBS eluates are already diluted to 1/50 after elution so only a 1/8 dilution is necessary

NB - We advise depositing 15µL of eluate in 105µL of dilution buffer.

For rodent dilute the DBS eluates to 1/400 in dilution buffer. The rodent DBS eluates are already diluted to 1/50 after elution so only a 1/8 dilution is necessary

NB - We advise depositing 15µL of eluate in 105µL of dilution buffer.

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For example, if 22 different recombinant proteins have been successfully coupled, the volume of the bead mix for one well will be $22 \times 0.25 = 5.5 \mu\text{L}$ for one well, which will be completed with 44.5 µL of dilution buffer.

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For 100 wells, 25 µL of each paired bead must be deposited per well, all in a final volume of 5000 µl qsp of dilution buffer.

Thus, if 22 different recombinant proteins have been successfully coupled, the volume of the mix of beads for a 96-well plate will be $22 \times 25 \mu\text{L} = 550 \mu\text{L}$ which will be completed with 4450 µL of dilution buffer.

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- Identify and mark the plate.
- Vortex the Mix of Marbles for 30 sec and pour it into a tank (provided for this purpose only)
- Add 50µl of bead mix/well to all wells with a multichannel pipette

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- Distribute 100 µl/well of diluted samples by mixing by aspiration/discharge with the multichannel pipette or 100 µl of buffer for the blank well.

- Replace the plastic lid on the plate and cover it with aluminum foil. Place it on the stirrer which has been placed at 4°C.
- Set the agitator to 400rpm.

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- Prepare the final conjugate solution in dilution buffer.
- Remove the stirrer plate and remove the aluminum (place the stirrer at room temperature)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 150µL/well of dilution buffer (Washing 1/3)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (2/3 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (3/3 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Off the support, deposit 50µL/well of the final conjugate solution
- Replace the plastic lid on the plate and place it on the stirrer and cover it with aluminum foil.
- Holding the plate firmly, set the stirrer at 1200 RPM for 1 minute to resuspend the beads.
- Reduce the stirrer speed to 400 RPM and incubate for **30 minutes at room temperature.**

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- **MagPix Calibration should be performed weekly when the MagPix is in use.**

- In the AutoMaint tab select Calibration Verification
- Fill the tank with 70% Isopropanol in the locations indicated on the screen by the AF zone.
- Vortex for 30s the 3 vials of beads from the verification kit and the vial from the calibration kit
- Add 5 to 6 drops of beads of each color in the wells of the bar indicated on the screen
- Open the hatch from the software using the "eject/retract" button
- Place the tank and the bar and click again on "eject/retract" to close the hatch.
- Click on "Start" to launch the "Calibration Verification" which lasts about 12 minutes.

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- **The MagPix Performance verification should be performed daily when the Magpix is in use**

- In the AutoMaint tab select the Performance Verification
- Fill the tank with 70% Isopropanol in the location indicated on the screen by the AF zone.
- Vortex for 30s the 3 vials of beads from the verification kit
- Add 5 to 6 drops of beads of each color in the wells of the bar indicated on the screen
- Open the hatch from the software using the "eject/retract" button
- Place the tank and the bar and click again on "eject/retract" to close the hatch.
- Click on "Start" to launch the "Performance Verification" which lasts about 5 minutes.

- **batch preparation (with an existing protocol)**

- in the batches tab click on "new batch from an existing protocol" and select "PanCoV protocol".
- Set the plate format (position of blank and samples)
- Import the plate plan (linear .txt format) and complete the dilution column
- In the "post batch routine" pane, select the "post batch routine" command.
- Fill the tank with distilled water, 70% Isopropanol and 0.1N soda (NaOH) in the locations indicated on the screen.
- Open the hatch from the software using the "eject/retract" button
- Place the tank and the bar and click again on "eject/retract" to close the hatch.
- Save the batch which then appears as pending (Pending)

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Streptavidin-PE (1mg/mL) is used at **1µg/mL final** (for a final volume of 50µL of solution/well).

Example for a plate: you need 100 wells x 50µL of solution at 1µg/ml = 5 ml of conjugate solution at 1µg/ml final (i.e. 5 µL of Streptavidin-PE at 1mg/ml) in 4995µL of dilution buffer)

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- Prepare the Streptavidin-PE solution at 1 µg/mL final in dilution buffer
- Cover the tube of diluted Streptavidin-PE with aluminum to protect it from light.
- Remove the agitator plate and remove the aluminum
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 150µL/well of dilution buffer (Washing 1/3)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (2/3 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Outside the support, add 100µL/well of dilution buffer (3/3 wash)
- Place the plate on the magnetic support for 2 min and wash by inversion
- Off the support, deposit 50µL/well of the Streptavidin-PE solution at 1 µg/mL final
- Replace the plastic lid on the plate and place it on the stirrer and cover it with aluminum foil.
- Holding the plate firmly, set the stirrer at 1200 RPM for 1 minute to resuspend the beads.
- Reduce the stirrer speed to 400 RPM and incubate for 10 minutes at room temperature.

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- Remove the agitator plate and remove the aluminum
- Place the plate on the magnetic support for 2 min and wash by inversion
- Out of support, deposit 150µL/well of **reading buffer (Washing 1/3)**
- Place the plate on the magnetic support for 2 min and wash by inversion
- Out of support, deposit 100µL/well of **reading buffer (Washing 2/3)**
- Place the plate on the magnetic support for 2 min and wash by inversion
- Out of support, deposit 100µL/well of **reading buffer (Washing 3/3)**
- Place the plate on the magnetic support for 2 min and wash by inversion
- Excluding support, deposit 50µL/well of reading buffer
- Replace the plastic lid on the plate and place it on the stirrer and cover it with aluminum foil.
- Holding the plate firmly, set the stirrer at 1200 RPM for 1 minute to resuspend the beads.
- Add 50µl/well of reading buffer to obtain a final volume of 100µl/well.
- Open the hatch from the software using the "eject/retract" button
- Remove the cover and place the plate on the MagPix support and click again on "eject/retract" to close the hatch
- Open the batch selecting it and clicking on Edit
- Click on the "START" button, and validate so that the reading of the plate starts

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- Fill the tank with mQ water and bleach diluted at 10% in the places indicated on the screen.
- Open the hatch from the software using the "eject/retract" button
- Place the tank and the bar and click again on "eject/retract" to close the hatch.
- Click on "Start" to launch the shutdown which lasts about 5 minutes.

- Turn off the device via the software then de-identify
- Close the Xponent software by clicking on Exit
- Set the switch on the back of the device to off
- Turn off the computer